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## INSTRUCTIONS FOR USE

Non-precious, beryllium and nickel free CoCr based bonding alloys for ceramic application regarding EN ISO 22674, Typ 5

### Nominal Alloy composition:

Co 63,5 | Cr 29,2% | Mo 5,0% | Si 1% | Mn, C < 1%

### Nominal Properties:

colour	silver
Melting interval	1320-1350 (°C)
Casting temperature	ca. 1450 (°C)
Vickers hardness	350 (HV10)
Elongation limit Rp 0,2	ca.215.000
Tensile strength Rm	600 (Mpa)
Modulus of elasticity	880 (Mpa)
Tensile elongation A5	6,5 (%)
Density	8,2 g/cm <sup>3</sup>

### Indications:

For the construction of frameworks according to the law for medical products. DIN EN ISO 22674 and rhe Directive 93/42/EEC on medical devices.

### Investing and Melting:

Use phosphate-bonded partial-denture investment materials  
 Use ceramic crucibles. Once melted, cast the molten alloy as quickly as possible. When melting with a flame, do not overheat the molten alloy. Do not use flux powder. Warming the ceramic crucible in advance simplifies casting. Use a separate crucible per alloy.

### Devesting:

Do not hit the button. Use 250 µm aluminium oxide for blasting off coarse remains of invest- ment material and 110 µm aluminium oxide for removing fine remains of investment.

### Trimming:

Never trim metal frameworks with diamond burs. When using abrasive stones, ensure that they are ceramic-bonded. The best results are achieved with tungsten carbide cutters.

### Soldering and welding:

With CoCr soldering rods Applicable using methods for non precious dental alloys and such as laser welding.

### Warranty:

Our recommendations for users are based on so-called approximate values determined in our test laboratories. These values can only be guaranteed if the procedures quoted are adhered to. The user bears the responsibility for processing the products. We accept no liability for defective results.

realloy

realloy e.K  
 Siemensdyk 50  
 D-47807 Krefeld  
 Germany

Website [www.realloy.net](http://www.realloy.net)  
 e-mail [info@realloy.net](mailto:info@realloy.net)  
 Tel + 49 (0) 2151 - 48 64 978  
 Fax + 49 (0) 2151 - 48 64 981