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INSTRUCTIONS FOR USE

Non-precious, beryllium and nickel free CoCr based bonding alloys for ceramic application regarding EN ISO:22674, Typ 5

Nominal Alloy composition:

Co 59,8-| Cr 24,9 % | W 6 % | Mo 4,8 % | Ga 3,0 % | Si 1,0 | Nb, Fe, Mn < 1 %

Nominal Properties :

| | |
|---------------------------|-----------------------|
| colour | white |
| Melting interval | 1310 - 1405 (°C) |
| Casting temperature | 1550 (°C) |
| Vickers hardness | 280 (HV10) |
| Modulus of elasticity | 225 (GPa) |
| Elongation limit (Rp 0.2) | 530 (Mpa) |
| Tensile strength (Rm) | 700 (Mpa) |
| Tensile elongation(A5) | 14 (%) |
| Density | 8.6 g/cm ³ |
| WAK (25 - 500 °C) | 14,0 (µm/mK) |
| WAK (25 - 600 °C) | 14,2 (µm/mK) |

Indications:

For fabricating dental crowns, bridges and frameworks according to the Directive 93/42/EEC on medical devices.

Waxing-up:

As non-precious alloys undergo high shrinkage after casting, to achieve optimum fit an approximately 0.1 mm thick spacer should be created using die spacer. The same effect can be achieved with vacuum-formed copings and spacer foil. For facing with porcelain, the metal must be at least 0.3 mm thick. For facing with composite and retention beads, the metal must also be at least 0.3 mm thick. Avoid sharp edges throughout the entire wax pattern.

Spruing:

Sprues attached directly to crowns should be 3.5-4.0 mm thick.

T-bar:

The connector between the T-bar and pattern should be 2.0-2.5mm thick. The connector between the crucible former and T-bar should be 3.5-4.0 mm thick. The junctures should be rounded and smoothed.

Investing:

Use phosphate-bonded investments and adhere closely to their manufacturers' instructions.

Melting:

Use ceramic crucibles. Once melted, cast the molten alloy as quickly as possible. When melting with a flame, do not overheat the molten alloy. Do not use flux powder. Warming the ceramic crucible in advance simplifies casting. Use a separate crucible per alloy.

Devesting:

Do not hit the button. Use 250 µm aluminium oxide for blasting off coarse remains of investment material and 110 µm aluminium oxide for removing fine remains of investment. The interiors of crowns should be sandblasted at a maximum pressure of 2 bars.

Trimming:

Never trim metal frameworks with diamond burs. When using abrasive stones, ensure that they are ceramic-bonded. The best results are achieved with tungsten carbide cutters.

Applying the porcelain:

All surfaces to be built up with porcelain must be sandblasted thoroughly. The oxide firing should be carried out in vacuum. After firing, sandblast again with 110 µm grit at a pressure of 2 bars. Then clean with steam, running water or an ultrasonic cleaner. Do not touch the surface with your fingers. Allow to cool normally after every firing. When building up the porcelain, adhere strictly to the porcelain manufacturers' instructions.

Soldering and welding:

With CoCr soldering rods Applicable using methods for non precious dental alloys and such as laser welding.

Warranty:

Our recommendations for users are based on so-called approximate values determined in our test laboratories. These values can only be guaranteed if the procedures quoted are adhered to. The user bears the responsibility for processing the products. We accept no liability for defective results.

realloy

realloy e.K
Siemensdyk 50
D-47807 Krefeld
Germany

Website www.realloy.net
e-mail info@realloy.net
Tel + 49 (0) 2151 - 48 64 978
Fax + 49 (0) 2151 - 48 64 981